


# ED COIL

Electro Detach Coil




# ED COIL 10 ExtraSoft *Optical*




More flexible and foldable than ever before.

## Feature

- The advanced flexibility of the coil realizes the Loop-in-Loop coil-folding feature which enables coil-packing without remaining open space.
- Wavy stretch resistant wire enables smooth coil-folding until the end.



2-4 ES (ø 2mm×4cm) Type RS tightly packed in a 2mm aneurysmal model. The highly flexible coil forms a small loop inside an outer loop (Loop-in-loop coil).



The stretch resistant structure created with 10µm platinum wave-processed stranded wire enables the coil to be delivered into small open space and be folded until the end.

Flexibility increased by 20% comparing to the non-Type RS ExtraSoft. The improved foldability and low shape memory enables the coil to access and fill the open space.


# ED COIL ∞ 10 [Infini]



A low shape memory coil provides flexible filling regardless of the size and shape of an aneurysm.

## Feature

- Conforming outward to reduce compartment space
- Compactly folding and entering into even a open space
- Conforms in accordance with the shape of an aneurysm even with irregular shape
- Provides a tight embolization with a small number of coils for a vessel occlusion



Infini Soft's large coil-diameter reduces the risk of coil-protrusion through a stent strut into the parent artery, and it enables its use for both framing and filling.

Waving Inner Wire Technology      Proximal Wider Pitch Technology



Waving Inner Wire Technology enables coil-folding in smaller size despite its large primary coil-diameter.



Infini Soft conforms outward to the aneurysm wall, thereby contributing to uniform packing even in an aneurysm with irregular shape.


# ED COIL 14 Standard *Optical*



αSpiral coil spreading randomly

## Feature

- The coil with 0.014 inch primary coil-diameter may provide sufficient coil embolization with a small number of coils
- The unique αSpiral coil structure provides random coil movement by deflection according to the aneurysm wall
- Strong framing by the unique αSpiral structure
- Suitable for both framing and filling
  - ≤ 6 mm (small pitch): suitable for filling
  - ≥ 7 mm (large pitch): suitable for both framing and filling



The random coil spreading with the αSpiral structure enables to change direction of the coil expanding in accordance with the shape of the wall of the aneurysm.



The unique αSpiral shape, designed to form a strong frame with random loops by having 1:1 ratio between coil diameter and pitch (≥ 7 mm coil)




αSpiral (φ10mm×L30cm) in a model of a 10 mm aneurysm




## The flexible pusher wire prevents the kick-back behavior of microcatheter and provides stable coil deployment.

- The tip (through 30mm from the distal end) is very flexible, therefore reducing the risk of kick-back of the microcatheter from the aneurysm.
- The coil is designed to minimize joints and attachments length, allowing the coil to be delivered even through the tortuous blood vessels.
- Stiffness of the pusher near the proximal part has doubled, allowing smoother delivery and achieving better response to the coil.



Easier to maintain catheter angle compared to other competitive products, allowing stable coil deployment even when the catheter tip is unstable in the aneurysm



Stiffness of the pusher near the proximal part is designed to optimize pushability and coil response.

## ELECTRO DETACH GENERATOR v4 *Detach Point Detection System*

The Electro Detach Generator v4 detects and indicates by sound and light signals the best Detach Point even when the position of the second marker of a microcatheter is not detectable under a X-ray fluoroscopy, Easy to use, disposable, sterilized, detachment of the coil can be done in single step operation by pressing 2 buttons simultaneously.

### Red lamp

The red lamp is lit while the detaching part of the ED COIL stays within the microcatheter.

### Green lamp

The green lamp is lit when the detaching part of the ED COIL has come out of the microcatheter.

### Orange lamp

The orange lamp is lit when a current leak or other trouble occurs.



### DETACH button

Detaching output is generated when the left and right buttons are pressed simultaneously while the green lamp is lit.

### Manual Detachment

When it is necessary to detach the ED COIL while the red lamp is lit, press and hold the left and right buttons simultaneously for longer than 2 seconds, and a detaching output is generated.

### RESET button

If the green lamp is lit while the detaching part is still within the microcatheter, press the reset button to restart detecting the detach point.

### POWER ON ribbon

Pull off the POWER ON ribbon and the Electro Detach Generator v4 is powered on.



The red lamp is lit when the ED clip is connected to the pusher.




The green lamp is lit when the detaching part of the ED COIL has come out of the microcatheter. Subsequently, detachment of the coil can be done by pressing 2 buttons simultaneously.



Patient-side clip to be used with Hypodermic Needle made of stainless-steel without resin-coating (20G to 22G)



# SET UP



## Product Specifications

ED COIL **Type RS**

Pusher Length 187cm

### ED COIL 10 ExtraSoft *Optical*

Catalogue Number	Primary Coil Diameter (inch)	Secondary Coil Diameter (mm)	Primary Coil Length (cm)
353-01H1RS	<b>0.010</b>	1.5	1
353-01H2RS			2
353-01H3RS			3
353-0202RS		2	2
353-0203RS			3
353-0204RS			4
353-0206RS			6
353-0208RS			8
353-02H3RS		2.5	3
353-02H4RS			4
353-02H6RS			6
353-0303RS		3	3
353-0304RS			4
353-0306RS			6
353-03H6RS		3.5	6
353-03H8RS			8
353-0406RS		4	6
353-0408RS			8

### ED COIL $\infty$ 10 ExtraSoft [Infini]

Catalogue Number	Primary Coil Diameter (inch)	Secondary Coil Diameter (mm)	Primary Coil Length (cm)
353-1610RS	<b>0.010</b>	16	10
353-1615RS			15

### ED COIL $\infty$ 10 Soft [Infini]


Catalogue Number	Primary Coil Diameter (inch)	Secondary Coil Diameter (mm)	Primary Coil Length (cm)
352-1620RS	<b>0.010</b>	16	20
352-1630RS			30

### ED COIL 14 Standard *Optical*

Catalogue Number	Primary Coil Diameter (inch)	Secondary Coil Diameter (mm)	Primary Coil Length (cm)	
350-0209RS	<b>0.014</b>	2	9	
350-0312RS		3	12	
350-0412RS		4	12	
350-0515RS		5	15	
350-0620RS		6	20	
350-0720RS		7	7	20
350-0730RS			30	
350-0820RS		8	8	20
350-0830RS			30	
350-0930RS			9	30
350-1030RS			10	30
350-1230RS		12	30	

## EDG

Catalogue Number
347-404



Do not re-use



Sterilized using ethylene oxide



Consult instructions for use

importato in esclusiva da:

Manufacturer

**KANEKA CORPORATION**

3-18, 2-Chome, Nakanoshima, Kita-ku, Osaka-city, OSAKA, 530-8288 JAPAN

Tel.No. : (+81)-(0)6-6226-5256 **CONTACT** Tel.No. : (+81)-(0)3-5574-8136

EC-Representative

**KANEKA PHARMA EUROPE N.V.**

Nijverheidsstraat 16, 2260 Westerlo-Oevel, BELGIUM

Tel.No. : (+32)-(0)14-256-297 Fax No. : (+32)-(0)14-256-298

**EUKON**  
Future Healthcare Now



**KANEKA  
MEDICAL  
PRODUCTS**