



medical **ECONET**
GERMANY

PROVIEW 10 / 12

Accessory list



Table of Contents

- Standard Accessory 2
- ECG 3
- SpO2 4
- Non-Invasive Blood Pressure 5
- Temperature 6
- Invasive Blood Pressure 6
- C.O. 7
- ETCO2 8
 - MicroFlow Stream (BLT)
 - Main Stream (BLT)
- Others 11

Standard accessory



5-lead Snap Connector ECG Cable

41.20-5362

IEC, 12 pin, AQS, Encrypted
Compatible with PROVIEW10, 12



Adult Disposable Adhesive Electrodes

70.10-3700

Snap Connector, 30pcs / pack
Compatible with PROVIEW10, 12



Analog SpO2 adult soft reusable sensor

41.20-5363

12 pin, Encrypted, One piece
Compatible with PROVIEW10, 12



NIBP PVC hose

41.20-5369

3 m, A series connector
Compatible with PROVIEW10, 12



Adult cuff

21.10-7051

with connector M5125, 250mm ~ 350 mm
Compatible with PROVIEW10, 12



Lithium Battery

41.20-5375

11.1 V / 2500mAh
Compatible with PROVIEW10, 12



User's Manual

Manual

220 V
Compatible with PROVIEW10, 12

ECG



5-lead Snap Connector ECG Cable

41.20-5362

IEC, 12pin, AQS, Encrypted
Compatible with PROVIEW10, 12



3-lead Snap Connector ECG Cable

41.20-5361

IEC, 12pin, AQS, Encrypted
Compatible with PROVIEW10, 12



Adult Disposable Adhesive Electrodes

70.10-3700

Snap Connector, 30 pcs/pack
Compatible with PROVIEW10, 12



Pediatric / Neonatal Disposable Adhesive Electrodes

70.10-2107

Snap Connector, 30 pcs/pack
Compatible with PROVIEW10, 12

SpO2

	Adult Reusable Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5363
	Adult Soft Reusable Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5366
	Extension Cable 9pin Compatible with PROVIEW10, 12	41.20-5368
	Pediatric Reusable Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5364
	Pediatric Soft Reusable Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5367
	Neonatal Reusable Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5365
	Adult / Neonatal Disposable Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5440
	Reusable Y Multi-Site Sensor 9pin, Encrypted Compatible with PROVIEW10, 12	41.20-5391

NIBP



NIBP PVC Hose

41.20-5369

3m

Compatible with PROVIEW10, 12



Adult cuff

21.10-7051

with connector M5124, 250 mm ~ 350 mm

Compatible with PROVIEW10, 12



Large Adult cuff

21.10-7052

with connector M5125, 330 mm ~ 470 mm

Compatible with PROVIEW10, 12



Pediatric Cuff

21.10-7050

with connector M5123, 180 mm ~ 260 mm

Compatible with PROVIEW10, 12



Neonatal Cuff

21.10-7037

with connector M5121, 60 mm ~ 110 mm

Compatible with PROVIEW10, 12

TEMP



Skin Temperature Probe

41.20-5371

for adult, pediatric, neonatal
Compatible with PROVIEW10, 12



Rectal Temperature Probe

41.20-5372

for adult, pediatric, neonatal
Compatible with PROVIEW10, 12

IBP



IBP Kit

41.20-5389


Disposable pressure transducer kit
(IBP cable 1pc, IBP disposable transducer 25 pcs)
Compatible with PROVIEW 12



IBP Cable

41.20-5388

with connector
Compatible with PROVIEW 12



IBP Disposable Transducer

21.11-5951

PT-01
Compatible with PROVIEW 12

C.O.



C.O. Module

41.20-5381

without Accessory Kit

Compatible with PROVIEW 12



C.O. Kits

41.20-5382

1 Cardiac Output cable
1 In-line Injection temperature probe (SP4042)
1 In-line Injection temperature probe housing (SP5045)
1 Control Syringe

Compatible with PROVIEW 12



Cardiac Output cable

41.20-5441

with connector

Compatible with PROVIEW 12



In-line Injection temperature probe

41.20-5442

SP4042

Compatible with PROVIEW 12



In-line Injection temperature probe housing

41.20-5443

BD SP5045

Compatible with PROVIEW 12



Control Syringe

41.20-5444

Compatible with PROVIEW10, 12

ETCO2

MicroFlow Stream



Capno_S ETCO2

41.20-5385


Water filter 2 pcs, Disposable CO2 sample tube 1pc, Airway adapter 1pc
Compatible with PROVIEW 12



12 Pin ETCO2 Module Extension Cable

41.20-5390

for Masimo Gas Module
Compatible with PROVIEW 12



Water Filter

41.20-5445

T4F
Compatible with PROVIEW 12



L type 3 way airway adapter

41.20-5446

for intubated patient
Compatible with PROVIEW 12



Disposable CO2 Nasal Cannula

41.20-5447

for Adult
Compatible with PROVIEW 12

ETCO2

Main Stream



Capno_M ETCO2

41.20-5386

with Disposable airway adapter for Adult / Pediatric 1pc
Compatible with PROVIEW 12



12Pin ETCO2 module extension cable

41.20-5390

for Masimo Gas module
Compatible with PROVIEW 12



Disposable Airway Adapter

41.20-5448

for Adult / Pediatric
Compatible with PROVIEW 12

Others



DM Module

41.20-5373

12pin, Encrypted
Compatible with PROVIEW10, 12



Mobile Stand

41.21-5374

Quick lock, RS002-5-BLT-S
Compatible with PROVIEW10, 12



Wall Mount

41.20-5450

Quick lock, WM001C-BLT-S
Compatible with PROVIEW10, 12



Lithium Battery

41.20-5375

11.1 V / 2500 mAh
Compatible with PROVIEW10, 12



Lithium Battery

41.20-5376

11.1 V / 5000 mAh
Compatible with PROVIEW10, 12



Printing Paper

41.20-5377

50 mm
Compatible with PROVIEW10, 12



WiFi Module

41.20-5378

Compatible with PROVIEW10, 12



Central Monitor

41.20-5392

Compatible with PROVIEW10, 12



medical ECONET GmbH
Im Erlengrund 20 46149 Oberhausen Germany
www.medical-econet.com

Distributor

*All specifications are subject to change without notice