


TECHNICAL SPECIFICATIONS

- PTCA Balloon Catheter Rapid Exchange System (RX)
- Non compliant balloon dilatation catheter
- Balloon material: Polyamide
- Balloon transition Angles: Tapered shoulders at 42°
- Balloon inflation pressures: Nominal 8 ATM/Bars.
Rated Burst Pressure (RBP): 20 ATM/Bars

- Catheter tip entry profile: 0.016"
- Catheter length: 140 cms
- Catheter shaft diameters: 1.9 FR (proximal) and 2.4 FR (distal)
- Compatible Guide wire: 0.014"
- Minimal Guiding catheter compatible: 5 FR (0.058mm)
- Kissing Balloon technique: 2 balloons on 6 FR guiding catheter




ORDERING INFORMATION

| BALLOON DIAMETER | BALLOON LENGTH | | | | | | |
|------------------|----------------|-----------|-----------|-----------|-----------|-----------|-----------|
| | 8 mm | 10 mm | 12 mm | 15 mm | 20 mm | 25 mm | 30 mm |
| 1.50 mm | NC-150-8 | NC-150-10 | NC-150-12 | NC-150-15 | NC-150-20 | - | - |
| 2.00 mm | NC-200-8 | NC-200-10 | NC-200-12 | NC-200-15 | NC-200-20 | NC-200-25 | - |
| 2.25 mm | NC-225-8 | NC-225-10 | NC-225-12 | NC-225-15 | NC-225-20 | NC-225-25 | NC-225-30 |
| 2.50 mm | NC-250-8 | NC-250-10 | NC-250-12 | NC-250-15 | NC-250-20 | NC-250-25 | NC-250-30 |
| 2.75 mm | NC-275-8 | NC-275-10 | NC-275-12 | NC-275-15 | NC-275-20 | NC-275-25 | NC-275-30 |
| 3.00 mm | NC-300-8 | NC-300-10 | NC-300-12 | NC-300-15 | NC-300-20 | NC-300-25 | NC-300-30 |
| 3.25 mm | NC-325-8 | NC-325-10 | NC-325-12 | NC-325-15 | NC-325-20 | NC-325-25 | NC-325-30 |
| 3.50 mm | NC-350-8 | NC-350-10 | NC-350-12 | NC-350-15 | NC-350-20 | NC-350-25 | NC-350-30 |
| 3.75 mm | NC-375-8 | NC-375-10 | NC-375-12 | NC-375-15 | NC-375-20 | NC-375-25 | NC-375-30 |
| 4.00 mm | NC-400-8 | NC-400-10 | NC-400-12 | NC-400-15 | NC-400-20 | NC-400-25 | NC-400-30 |
| 4.50 mm | NC-450-8 | NC-450-10 | NC-450-12 | NC-450-15 | NC-450-20 | - | - |
| 5.00 mm | NC-500-8 | NC-500-10 | NC-500-12 | NC-500-15 | NC-500-20 | - | - |

Catalogue Nr: 22-0040-8 Rev:00

NC CLEVER

NON COMPLIANT
PTCA BALLOON CATHETER



Conic Vascular Technology S.A.
Via Carlo Maderno, 23
6901 Lugano - Switzerland

Phone (+41) 91 910 77 88
Fax (+41) 91 923 91 19
info@conic.ch

www.conicvascular.com

Distributed by:

