


Ordering information

REF	Frechsize	Effective length (cm)	O.D.	I.D.	Guide segment (cm)	Side hole
20505603	5F	150	0.064"	0.056"	35	--
20606603	6F	150	0.078"	0.066"	35	--

Side-hole



REF	Frechsize	Effective length (cm)	O.D.	I.D.	Guide segment (cm)	Side hole
20505604	5F	150	0.064"	0.056"	35	4
20606604	6F	150	0.078"	0.066"	35	4



Please contact us for more information:

EUKON S.r.l.

Address: Via Claudio Coccia, 4 - Caserta CE - ITALY

Post code: 81100

E-mail: info@eukon.it

Tel: +3908231455909

Fax: +3908231455908

www.eukon.it


PN: ENG-Guiding Extension Catheter-210X285X4P-20180410




Expressman™ Guiding Extension Catheter

Embedded structure

Embedded structure designed provides safer operation support.



Gradually changed distal O.D. for extending your reaching.



4 side-hole, 90° spatial distribution

Auxiliary pathway while treating complex or long lesion.

