


Ordering information

REF	Frechsize	Effective length (cm)	O.D.	I.D.	Guide segment (cm)	Side hole
20505603	5F	150	0.064"	0.056"	35	--
20606603	6F	150	0.078"	0.066"	35	--

Side-hole



REF	Frechsize	Effective length (cm)	O.D.	I.D.	Guide segment (cm)	Side hole
20505604	5F	150	0.064"	0.056"	35	4
20606604	6F	150	0.078"	0.066"	35	4

CE₀₁₂₀

Please contact us for more information:


Address:
 Post code: E-mail:
 Tel: Fax:
 PN: ENG-Guiding Extension Catheter-210X285X4P-20180410




Expressman™
Guiding Extension Catheter

Embedded structure

Embedded structure designed provides safer operation support.



Gradually changed distal O.D. for extending your reaching.



4 side-hole, 90° spatial distribution

Auxiliary pathway while treating complex or long lesion.

