



TrueCross™

Single-use Microcatheter

*Redefine flexibility
and maneuverability*



Tapered Design Distal Shaft and Hydrophilic Coating


- Enhanced crossability and pushability.
- Low profile, easier to cross tortuous lesions and fine vessels, even in complex PCI and CTO procedures.

DLF (Double Layer Fusing)™ Technology

- Tapered and radiopaque soft tip, can be easily located during the operations.
- Excellent crossability in the lesions and tracking ability for the guide wires.

HybriBraiding™ Technology

- Unique body shaft design with high flexibility and torqueability.
- Maintaining lumen integrity in tortuous anatomy.



TrueCross™ Single-use Microcatheter Series

	TrueCross™ Pro		TrueCross™ Duo		
Model	03-19F135-01	03-19F150-01	03-21F135-01	03-21F150-01	03-21F165-01
Proximal Outer Diameter	0.87mm(2.6F)	0.87mm(2.6F)	0.91mm(2.7F)	0.91mm(2.7F)	0.91mm(2.7F)
Distal Diameter	0.64mm(1.9F)	0.64mm(1.9F)	0.71mm(2.1F)	0.71mm(2.1F)	0.71mm(2.1F)
Working Length	135cm	150cm	135cm	150cm	165cm

For more information, please visit our website at www.insight-med.com or contact us:



 **Insight Lifetech Co., Ltd.**
 Zone E, 3rd Floor, Building 3, Tinwe Business Park,
 No.6 Liufang Road, Baoan District, 518000
 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
 Tel: +86 (755) 23074845
 Email: export@insight-med.com



TrueCross is trademark of Insight Lifetech Co., Ltd. All data and photos on file.
 Illustrations are not to be considered as engineering drawings of photos.
 ©2023 Insight Lifetech Co., Ltd.